

Supplementary Materials: Stencil Printing—A Novel Manufacturing Platform for Orodispersible Discs

Henrika Wickström *, Rajesh Koppolu, Ermei Mäkilä, Martti Toivakka and Niklas Sandler

Table S1. Repeatability, reproducibility and sample stability of the HPLC method.

Repeatability (<i>n</i> = 6 within 1 day)			
Conc. ¹	10	30	50
µg/ml ²	10.11	32.12	49.64
Std	0.06	0.12	0.18
RSD%	0.56	0.41	0.35
Reproducibility (<i>n</i> = 12 within 3 days)			
Conc. ¹	10	30	50
µg/ml ²	10.13	29.92	49.93
Std	0.02	0.18	0.25
RSD%	0.22	0.58	0.5
Sample stability (<i>n</i> = 6 within 2 days)			
Conc. ¹	10	30	50
µg/ml ²	10.07	30.01	49.66
Std	0.05	0.15	0.01
RSD%	0.47	0.49	0.03

¹ theoretical and ² measured concentrations.