

Supplementary Materials

Can Aluminum Tolerant Wheat Cultivar Perform Better under Phosphate Deficient Conditions?

Table S1: List of primers use for quantitative RT-PCR (qRT-PCR).

Genes	Forward primer sequence (5'–3')	Reverse primer sequence (5'–3')	Accession number
<i>TaIPS1</i>	GACACTGAAGACTCGCACCA	CGGCGACTTCTCACCTCTAC	EU753151.1
<i>TaSOD2</i>	TGCCAGTGGAGATGTGTTTGT GAT	AAGCTGGTCTTTCCTTCCTGAT	AK333217.1
<i>TaSPX3</i>	GTGGAAGGACGAGTTCCTGAG C	TCCCGGTGTGTGATGATGAAGAA	AK333013.1
<i>TaAPT1</i> (Reference gene)	CGAATCAGTATGAAACAAGTT GTGACTCTT	TCCCAAACAGTTCAGAAAGGACAAAT	U22442.1